Subject: Deleting 3D Confomers

Posted by umaed on Fri, 14 Jul 2023 10:53:07 GMT

View Forum Message <> Reply to Message

Hello

I am trying to delete excess confomers I have made. These clutter the detail area. I add new chemicals to my database, but then want to generate new confomers to make sure the whole database has them. Is there a way to delete old confomers across the database to only have the new ones? (see pic)

Can I do confomers for one molecule at a time if I add one or 2 to the database?

Many thanks

File Attachments

1) Screenshot 2023-07-14 114712.png, downloaded 264 times

Subject: Re: Deleting 3D Confomers

Posted by thomas on Fri, 14 Jul 2023 10:59:11 GMT

View Forum Message <> Reply to Message

Data->Delete Columns...

Then select the 3D-coordinate columns you want to delete and press OK

Subject: Re: Deleting 3D Confomers

Posted by umaed on Fri, 14 Jul 2023 11:01:11 GMT

View Forum Message <> Reply to Message

Amazing! Seems obvious now!

Thank you